Resin-based Dental Restorative Material

ESTELITE Σ QUICK SYRINGE**ESTELITE Σ QUICK** PLT

CE 1639

UK 0120
CA

Instructions for Use

Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com/>

ENGLISH	2	БЪЛГАРСКИ	7	FRANÇAIS	14
ΕΛΛΗΝΙΚΑ	20	ITALIANO	27	PORTUGUÊS	33
ROMÂNĂ	39	ESPAÑOL	45	TÜRKÇE	52
HRVATSKI	58	SLOVENŠČINA	63		

This product is certified as a medical device under the legislation of the European Union on medical devices by SGS CE1639, exclusively for the indication(s) shown in this IFU. Other non-medical uses ascribed to this product are not within the scope of CE certification, and users should be aware product performance and/or safety has not been evaluated by SGS for those purposes.

Read all information, precautions and notes before using.

■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

- 1) ESTELITE Σ QUICK is a light-cured, radiopaque, submicron filled composite resin for use in anterior and posterior restorations and is indicated for all carious classes. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. ESTELITE Σ QUICK contains 82% by weight (71% by volume) of silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. All inorganic filler contained in ESTELITE Σ QUICK is a spherical submicron filler (mean particle size: 0.2 μm, particle size range: 0.1 to 0.3μm) that is designed to deliver excellent gloss retention and wear resistance. The monomer matrix contains Bis-GMA and Triethylene glycol dimethacrylate.
- 2) ESTELITE Σ QUICK implements Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE Σ QUICK can be shortened by one-third using a light curing unit with a wavelength range of 400 to 500nm (peak: 470nm) compared to our conventional resin composite. Please see the table depicting the relationship between curing time and increment depth (■ INDICATIONS FOR FILLING AND CURING).
- 3) RAP technology facilitates a longer working time for ESTELITE Σ QUICK.
- 4) ESTELITE Σ QUICK is provided either in SYRINGE or Pre Loaded Tip(PLT).

■ SHADE

- A1, A2, A3, A3.5, A4, A5, B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, OPA2, BW (Bleach White), WE (White Enamel) and CE (Clear Enamel/Incisal).
- Opalescent shades (OA1, OA2, OA3) have adequate opacity designed for blocking out the dark shine that comes through the oral cavity (Class III and Class IV restoration). The opalescent shades also serve as a dentin layer for the multi-shade layering technique. Shades OA1, OA2 and BW closely match those of primary teeth. However, they are not designed for masking metal or a darkly stained area such as tetracycline stained.
 - OPA2 is designed for masking slight stain or reconstructing a highly opaque tooth (use as a thin layer at the lingual cavity wall of class III and IV).
 - Please check with your local distributor for available shades.

■ INDICATIONS

- Direct anterior and posterior restorations including the occlusal surface and all carious classes
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

■ CONTRAINDICATIONS

ESTELITE Σ QUICK contains methacrylic monomers. DO not use ESTELITE Σ QUICK for patients allergic to or hypersensitive to methacrylic and related monomers.

■ PRECAUTIONS

- 1) DO NOT use ESTELITE Σ QUICK for any purpose other than those listed in these instructions. Use ESTELITE Σ QUICK only as directed herein.
- 2) ESTELITE Σ QUICK is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- 3) DO NOT use ESTELITE Σ QUICK if the safety seals are broken or appear to have been tampered with.
- 4) If ESTELITE Σ QUICK causes an allergic reaction or oversensitivity, discontinue its use immediately.
- 5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE Σ QUICK to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE Σ QUICK comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 6) Avoid contact of ESTELITE Σ QUICK with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE Σ QUICK comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE Σ QUICK comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.
 - If ESTELITE Σ QUICK comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
- 7) ESTELITE Σ QUICK should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.

- 8) To avoid the unintentional ingestion of ESTELITE Σ QUICK, do not leave it unsupervised within the reach of patients and children.
- 9) Clean the placement instruments with alcohol after use.
- 10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

■ INDICATIONS FOR FILLING AND CURING

Be sure to light-cure ESTELITE Σ QUICK extraorally and check the time needed for complete hardening of ESTELITE Σ QUICK with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light type	Intensity (mW/cm ²)	Curing time (seconds)	Increment depth (mm) ¹⁾	
			A1, A2, A3, A3.5, B1, B2, B3, C1, C2, OA1, OA2, CE(Inc.), BW, WE	A4, A5, B4, C3, OA3, OPA2
Halogen	400	10	1.9	1.6
		20	2.2	1.9
	600	10	2.0	1.7
		20	2.5	2.0
	800	5	1.8	1.6
		10	2.2	1.9
LED	260	10	1.7	1.4
		20	2.0	1.7
	900	5	1.7	1.5
		10	2.0	1.8
Plasma Arc	950	3	1.9	1.5

- 1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 “depth of cure” of ISO4049.

■ SPECIAL NOTES FOR THE USE OF PLT

- 1) PLTs are designed for single patient use only. To avoid cross infection, Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
- 2) Dispensers are not provided for ESTELITE Σ QUICK. Use a dispenser that fits the PLT of

ESTELITE Σ QUICK. For proper use and proper disinfection, see manufacturer's instructions.

- 3) Apply the composite at normal room temperature (18 - 30°C / 62 - 84°F). Material that is too cold might be difficult to extrude.
- 4) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■ STORAGE

- 1) Store ESTELITE Σ QUICK at a temperature between 0 - 25°C (32 - 77°F).
- 2) AVOID direct exposure to light and heat.
- 3) DO NOT use ESTELITE Σ QUICK after the indicated date of expiration on the syringe or PLT package.

■ DISPOSAL

To safely dispose of excess ESTELITE Σ QUICK, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

Follow local instructions for disposal.

Refer to the packaging waste disposal guide shown on the end of this IFU.

■ CLINICAL PROCEDURE

1. Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Shade Selection

Select the appropriate shade of ESTELITE Σ QUICK prior to placing the rubber dam using the shade guide provided for ESTELITE Σ QUICK.

- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

3. Isolation

A rubber dam is the preferred method to isolate the tooth from contamination.

4. Cavity Preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- Scalloped (undulant) bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL BASED MATERIALS** to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE Σ QUICK.

6. Bonding system

Apply a bonding system according to manufacturer's instructions.

7. Dispensing

7-1. PLT

Please read ■ **SPECIAL NOTES FOR THE USE OF PLT** before use.

- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into cavity or indirectly on the mixing pad.

7-2. SYRINGE

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. Re-cap the syringe immediately after dispensing.

- **DO NOT** apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring

Fill the cavity incrementally. Increments should not exceed the indicated depth (please refer to the aforementioned table).

- **DO NOT** mix ESTELITE Σ QUICK with other brands of resin composite or other shades of ESTELITE Σ QUICK to avoid incomplete cure or voids.

9. Curing

Light-cure each increment for at least the indicated time (please see aforementioned table),

keeping the curing light tip within a distance of 2 mm from the increment.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGo™ (DENTSPLY/Caulk), Sof-Lex™ (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINE™ Hybrid Diamond (COSMEDENT, Inc.).

* PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

- Instruct the patient adequately on the maintenance of the restoration after treatment.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPORTANT NOTE: The manufacturer of ESTELITE Σ QUICK is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.

Product specifications of ESTELITE Σ QUICK are subjected to change without notification. When the product specification change, the instruction and precautions may change also.

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Packaging waste disposal guide

(According to European Commission Decision 97/129/EC)

ENGLISH

Packaging waste disposal guide

SYRINGE (cap and body)	Dispose of as industrial waste.
PLT (cap and body)	Dispose of as industrial waste.
Packaging bag for PLT ESTELITE Σ QUICK PLT	PP5 (Polypropylene): Recycling. Follow local instructions for disposal.
Packaging bag for PLT ESTELITE Σ QUICK PLT Mini	LDPE4 (Low density polyethylene): Recycling. Follow local instructions for disposal.
Paperbox and/or partation	PAP21 (Non-corrugated fibreboard): Recycling. Follow local instructions for disposal.
IFU and/or guide	PAP22 (Paper): Recycling. Follow local instructions for disposal.

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* PoGo, Sof-Lex, Identoflex in D-FINE niso registrirane blagovne znamke Tokuyama Dental Corporation.